

Labycare-SEED COUNTER COMPUTERISED (SKSS-SLE2024)



Seed Counter Computerised Manufacturer by LABYCARE

Description :

This machine is designed for the quick determination of seed quantity in terms of numbers. It utilizes an advanced microcomputer to control its operations. The machine offers additional features such as a removable pan, batch setting, count setting, repeat setting, and the ability to enter commodity names

<u>Specifications :</u>	
Capacity	Counts maximum up to 999999 seeds in a lot.
Range	0.3mm to 15mm.
Memory	To store settings of 99 commodity parameters which include seed size, seed name, lot no fast and slow vibration strength and counts.
Controller	Micro Controller based.
Keypad	Membrane type.
Display	16x2 line LCD display with backlit. Seed name, size, count and vibration strength displayed on LCD while counting.
Feeding	Electromagnetic vibratory feed.
Removable pan	Knob provided for removing pan easily. Separate pan for seed size below 1mm.
Speed control	Soft Touch speed control no knob required.
In flight error	Programmable stop to

	compensate in flight error.
Pan	Metallic pan of volume-200cc.
Interface	RS 232 and parallel port interface for attaching computer and printer.
Additional pan	Additional pan for seed size 0.3 to 0.4mm (optional).

labycare

www.labycare.com | Email: labycare@gmail.com