

Labycare-MARVIN SEED/GRAIN ANALYZER (SKSS-SLE2029)



MARVIN SEED/GRAIN ANALYZER Manufacturer by LABYCARE

Description :

The machine enables fast and accurate determination of Thousand Grain Weight (TGW) and high-speed seed counting with an error margin of less than 1%. It also performs 2D seed dimension analysis, distinguishing main seeds from off-type seeds and debris. With an analysis time of just 1-3 seconds, it ensures rapid and efficient operation. The double-check function enhances accuracy, while its software-based system consolidates all checks into a single operation. Reliable data storage allows for saving counting and measurement results, with direct data transfer to Win Access, Excel, and user databases for seamless record management

Specifications :

- Determination of the thousand grain weight (TGW).
- Counting of seeds with high speed and accuracy.
- Determination of 2D dimensions of the seeds.
- Distinction between main seeds and off type seeds, debris etc.
- Very fast analysis (1-3s) and high counting accuracy (error < 1%).
- Double check in a jiffy, improves the accuracy.
- Software based, all checks in one operation.
- Reliable determination, storing of counting and measuring results.
- Direct data transfer to Win Acces/excel, user databases possible.

labycare

www.labycare.com | Email: labycare@gmail.com