

Labycare-360° Grain Analyzer (SKSS-SLE2030)



360° Grain Analyzer Manufacturer by LABYCARE

Description :

The Grain Analyser 360° is an automated inspection system designed to enhance and standardize physical purity quality control in the agribusiness sector. Utilizing advanced computer vision technology, it performs comprehensive analyses of grain samples, providing objective and traceable results

<u>Specifications :</u>	360° Inspection	Captures three images per seed to analyze the entire
	3D Granulometry	Measures the length, width, and thickness of each gra
	Defect and Foreign Body Detection	Identifies and classifies defects and foreign materials v
	Automated Classification	Classifies each grain and estimates both the total weig
	Accurate Reporting	Generates comprehensive reports containing data, classification, and real-size images.
	Processing Speed	Analyzes a 500 g sample in approximately 60 seconds
	Database Storage	All obtained information, along with product images, is

labycare

www.labycare.com | Email: labycare@gmail.com