

## Labycare-DEHUSKED GRAIN TESTING KIT (SKSS-SLE2034)



Dehusked Grain Testing Kit Manufacturer by LABYCARE

### **Description :**

This kit enables the identification of dehusked rice grains using a double staining method. Grains are sequentially treated with methylene blue and Metanil yellow, producing a color distinction—green for dehusked grains and yellow for others—allowing accurate counting in a 5 g rice sample.

<b><u>Specifications :</u></b>	Metanil yellow	25gm.
	Methyle blue	25gm
	HCL	100 ml
	Regent bottle	500ml – 03 nos.
	Measuring cylinder glass	50cc, 20cc, 10cc-one each.
	Glass flask	1 ltr. – 01 nos.
	Petri dish	4” – 04nos.
	Forcep	01 nos.
	Pipette	10ml – 01nos.
	Magnifier	04” – 01 nos.

---

labycare

www.labycare.com | Email: labycare@gmail.com