

## Labycare-GRAIN ANALYSIS KIT ( SKSS-SLE2035)



Grain Analysis Kit Manufacturer by LABYCARE

### **Description :**

This kit is a compact tool frequently used during procurement to check the quality of food grains. It can also be used with a seed grader to assess the fractions of various impurities, such as shrunken and immature grains. All items are packed in a briefcase for convenience. The kit is applicable to all types of seed grains, including cereals, oilseeds, pulses, vegetables, and forestry seeds.

It is designed to measure the fractions of grain lots, including shrunken, shriveled, broken grains, mud balls, stones, infected grains, undersized seeds, light chaff, straw, immature grains, and hollow grains.

<b><u>Specifications :</u></b>	Mechanical weighing balance	Upto 100g
--------------------------------	-----------------------------	-----------

Weight Box	Upto 100g
Enameled Plate	9" dia
Parkhi	9" Length
Sieve	3.30 Mesh
Magnifying Glass	4"
Scoop	Stainless steel
Measuring Cylinder	20cc
Petri Dish	Good quality.
Forcep.	Good quality.
Bag	These items are packed in a specially designed canvas bag.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)