

## Labycare-SEED TESTING KIT ( SKSS-SLE2036)



seed testing kit Manufacturer by LABYCARE

### Description :

A Seed Testing Kit is a portable diagnostic tool designed to assess seed quality by evaluating germination rate, viability, moisture, purity, and physical attributes. Ideal for agricultural labs and field use, it ensures reliable quality control of seeds before planting or sale, supporting higher crop performance and compliance with standards.

### Specifications :

- Compact, portable and fully equipped kit.
- Precision Digital weighing Balance for grain sample analyzing .
- Purity work board with lens for physical purity testing.
- Sampling parkhi
- Petri Dishes
- Germination Towel paper for germination test.
- Wax and butter papers for germination test.
- Equipped in the compact wooden box.

---

labycare

www.labycare.com | Email: labycare@gmail.com