

# Labycare-Hectoliter Weight Apparatus (SKSS-SLE2037)



Hectoliter Weight Apparatus Manufacturer by LABYCARE

## **Description :**

A Hectoliter Weight Apparatus is a precision instrument used to measure the bulk density (hectoliter weight) of grains like wheat, rice, barley, and pulses. It determines the weight per 100 liters of grain, helping assess quality and grading for trade, storage, and processing in agricultural and grain testing facilities.

## **Specifications :**

- Smooth finish of inner portion of apparatus.
- Digital balance, capacity 2000g, readout 1g.
- Supplied with a precise 500 ml calibrated container featuring a slot for inserting a measuring plate to ensure uniform volume.
- Apparatus made of high-quality stainless steel.
- Provided with hectolitre conversion chart.
- Very user friendly in operation.

---

labycare

www.labycare.com | Email: labycare@gmail.com