

Labycare-Seed Riffle Type Divider (Soil Type) (SKSS-SLE2041)



Seed Riffle Type Divider (Soil Type) Manufacturer by LABYCARE

Description :

A Seed Riffle Type Divider (Soil Type) is a sampling tool used to split soil-attached seed or small grain samples into uniform subsamples. Its riffle channels ensure even distribution, aiding representative selection for quality tests like viability, purity, and moisture. Ideal for lab and field use in seed and soil analysis.

Specifications :

Model	No. of Chutes	Width of Chutes	Hopper Size(L x W)	Pan	Construction
RT-12	12	1.2 cm	28 x 20 cm	2 no.	Stainless Steel
RT-14	14	1.2 cm	30 x 20 cm	2 no.	Stainless Steel
RT-24	24	1.2 cm	42 x 20 cm	2 no.	Stainless Steel

labycare

www.labycare.com | Email: labycare@gmail.com