

## Labycare-Seed Hand Test Sieves (SKSS-SLE2042)



Seed Hand Test Sieves Manufacturer by LABYCARE

### **Description :**

Seed Hand Test Sieves are manual screening tools used to separate seeds and grains by size for quality evaluation. With standardized mesh openings, they help determine seed size distribution, remove debris, and aid purity analysis. Ideal for agricultural labs and field sampling, they support accurate seed grading and quality assessments.

### **Specifications :**

- Brass Top Lid & Bottom Pan for Dia. 12" x 3" Sieves of zinc plated MS, Brass frame.
- Brass Top Lid & Bottom Pan for Dia. 8" x 2" Sieves of zinc plated MS, Brass frame.
- Brass frame with lid & pan. Screen size 4, 3.35, 1.7, 1mm, pan dia 12" x 3".
- Brass frame with screen as per customer required perforation size 9" x 2".

---

labycare

www.labycare.com | Email: labycare@gmail.com