


Labycare-Hand Test Sieves (Wooden) (SKSS-SLE2043)



Hand Test Sieves (Wooden) Manufacturer by LABYCARE

Description :

Hand Test Sieves (Wooden) are manual screening tools with wooden frames and mesh screens used to separate and classify materials by particle size. Ideal for seeds, grains, and soil samples, they help determine size distribution and remove debris, aiding accurate quality evaluation in agricultural labs and field testing.

<u>Specifications :</u> WOODEN SQUARE	: Set of 5 Sieves with Bottom Pan, Top Pan and Tie Rods for Paddy Sizes : 2.75 OH, 2.0 OH, 1.85 OH, 1.75 OH, 2.25 OH. 6818
	: Tie Rod & Bottom Pan for 30 x 30 cm (Sieve).

labycare

www.labycare.com | Email: labycare@gmail.com