

## Labycare-Seed Sieve Shaker (SKSS-SLE2044)



Seed Sieve Shaker Manufacturer by LABYCARE

### **Description :**

A Seed Sieve Shaker is an automated laboratory instrument used to vibrate and agitate stacked sieves for efficient particle size separation of seeds, grains, and powders. It ensures fast, consistent, and reproducible sieving results, improving accuracy in quality analysis, grading, and research applications in seed testing laboratories and processing facilities.

<b>Specifications :</b> Type	: Batch type.
Vibration	: Adjustable.
Voltage	: 220V, single phase.
Sieve size	: 13"
Jally dia (Round Hole)	1.60mm RH, 2.4mm RH, 2.80mmRH, 3.5mmRH, 4.0mm RH, 6.0mm RH
Jally dia (Oblong Hole)	1.85mm OH, 2.25mm OH, 2.50mmOH , 2.

---

labycare

www.labycare.com | Email: labycare@gmail.com