

## Labycare-Computerised Hot Air Seed Dryer Cabinet Type (Single Chamber) (SKSS-SLE2046)



Computerised Hot Air Seed Dryer Cabinet Type (Single Chamber) Manufacturer by LABYCARE

### Description :

The cabinet features a triple-walled construction with gasketed doors. The inner chamber is made of stainless steel, while the outer layer is constructed from mild steel with a painted finish. The trays inside the chamber are adjustable. Tubular heaters are mounted between the middle sheet and the inner chamber on both sides. A motorized blower, designed for air circulation, is located below the chamber and is made of perforated stainless steel. The unit is mounted on a trolley for easy portability.

The doors are equipped with a high-quality locking system to ensure the security of the samples under test. Adjustable air ventilators are provided on both sides of the dryer. The ventilator holes can be adjusted to regulate the internal temperature by sliding the knobs up and down.

<b>Specifications :</b> Capacity	34 cu. ft.
inside Chamber	Stainless steel
Inside Chamber Size	90 cm x 90 cm x 120 cm (approx.) (L x B x H)
Outer chamber	Mild steel with paint finish
Electronics	PID Control, Microcontroller based Digital display of temperature.
Temperature	Ambient plus 5°C to +60°C ±2°C.
Alarm	Audiovisual alarm for programmable deviation.
Display	16x2 line alphanumeric LCD display of temperature
Computer Interface	RS 232 Serial Port.
No. of Trays	2
Tray Size (L x b x h)	90x89x5cm deep
Operating Voltage	220V ± 10% AC, 50Hz.

---

labycare

www.labycare.com | Email: labycare@gmail.com