

Labycare-Seed-Grain Sampling Trier (SKSS-SLE2047)



Seed-Grain Sampling Trier Manufacturer by LABYCARE

Description :

The Sampling Trier is designed to extract accurate samples from a given lot, ensuring maximum precision. It improves sample quality and prevents physical mixing of samples.

<u>Specifications :</u>	Type / Length	Outside Diameter	No. of Openings
	Double Tube, 75 cm	1.25 cm	5
	Double Tube, 90 cm	2.5 cm	5
	Double Tube, 150 cm	3.75 cm	5

labycare

www.labycare.com | Email: labycare@gmail.com