

## Labycare-Seed Blower - Dakota Type (SKSS-SLE2049)



Seed Blower - Dakota Type Manufacturer by LABYCARE

### **Description :**

Seed Blower (Dakota Type) is a laboratory instrument used to separate light impurities such as chaff, dust, and empty seeds from heavier seeds using controlled airflow. It helps in determining seed purity and quality during testing. The adjustable air control ensures accurate and efficient separation of seeds in laboratories.

<b><u>Specifications :</u></b> Function	Used for purity testing, docking & cleaning of grains
Timer facility	To set cleaning time.
Blower speed controller	Through mechanical valve.
Portable	Through caster wheels.
Transparent acrylic pipe	3 types diameter 1.5?, 3? & 4?
Motor	Single phase electric motor 0.5 HP.
Blower housing	Made of aluminium.

---

labycare

www.labycare.com | Email: labycare@gmail.com