

Labycare-Laboratory Aspirator - Bates Type (SKSS-SLE2050)



Laboratory Aspirator - Bates Type Manufacturer by LABYCARE

Description :

Laboratory Aspirator (Bates Type) is used in seed testing laboratories to separate light impurities such as chaff, dust, and broken particles from seed samples through controlled air suction. It ensures accurate purity analysis and grading of seeds. Adjustable airflow allows efficient separation while maintaining the integrity of good quality seeds.

<u>Specifications :</u>	Hopper capacity	2 kg approx.
	Motor power	1/10HP.
	Pan	Triangular for final collection
	Jar	Transparent for light impurities.
	Adjustment	Air & feed rate.

labycare

www.labycare.com | Email: labycare@gmail.com