

Labycare-Digital Temperature & Humidity Indicator (SKSS-SLE2055)



Digital Temperature & Humidity Indicator Manufacturer by LABYCARE

Description :

Digital Temperature & Humidity Indicator is used to measure and display temperature and relative humidity levels in laboratories, seed storage rooms, and controlled environments. It provides accurate digital readings for monitoring environmental conditions. The device helps maintain proper storage conditions to preserve seed quality and prevent moisture-related damage.

Specifications :

HUMIDITY

Range	5 to 95% RH
Resolution	0.1% RH
Accuracy	±2% RH at 25°C

TEMPERATURE

Range	-40°C to +60°C
Resolution	0.1°C
Accuracy	±0.5°C at 25°C
Power	230V ±10%, 50 Hz

labycare

www.labycare.com | Email: labycare@gmail.com