

Labycare-INDENTED CYLINDER (SKSS-SLE2059)



Indented Cylinder Manufacturer by LABYCARE

Description :

An Indented Cylinder is a laboratory and agricultural testing device used to separate and grade seeds or grains based on length. The rotating cylinder contains indents that pick up shorter grains while longer grains remain outside. It is widely used in seed testing laboratories for purity analysis and seed quality evaluation

<u>Specifications :</u>	Capacity (wheat basis)	1 Ton/hr	2Tons/hr	3Tons/hr	4
Motor (3 phase 440V)		01HP	2HP	2HP	3
Cut/Short Seeds Outlet		By auger	By auger	By auger	E

labycare

www.labycare.com | Email: labycare@gmail.com