

## Labycare-MIST-O-MATIC SEED TREATER (SKSS-SLE2062)



Mist-O-Matic Seed Treater Manufacturer by LABYCARE

### **Description :**

The Mist-O-Matic Seed Treater is a laboratory and seed processing device used to apply fungicides, insecticides, or protective chemicals uniformly on seeds. It uses a fine mist spray system to ensure even coating without damaging seeds. Widely used in seed testing laboratories and agricultural research for accurate seed treatment.

### **Specifications :**

- Seed Measuring and Feeding Unit:
- Hopper Capacity : 100 kg (depends on seed size and density)
- Level is controlled by a "High accurate precision Sensor".
- Machine housing with rubber or PU moulded Peddler is used for accurate feeding and to eliminate any possibility of seed damaged during feeding.
- High precision Geared Motor for Peddler movement.

---

labycare

www.labycare.com | Email: labycare@gmail.com