

Labycare-Kymograph Various, (Kymograph Acceseries) (SKSS-MD2001)



Kymograph Various, (Kymograph Acceseries) Manufacturer by LABYCARE

Description :

Kymograph is a laboratory instrument used to record physiological movements such as muscle contractions, heartbeat, and respiration. It consists of a rotating drum with recording paper and stylus arrangement. Widely used in educational and research laboratories, it ensures accurate graphical recording and comes with essential accessories for versatile experiments.

Specifications :

Available in manual (spring-driven) and electric motor models
Durable metal/aluminum rotating drum (approx. 150–200 mm length)
Multi-speed control with gear or stepper mechanism
Stylus-based recording on smoked or plain paper
Stable heavy cast iron base for vibration-free operation
Electric model operates on 220–240V AC supply

labycare

www.labycare.com | Email: labycare@gmail.com