

Labycare-Spirometer Various (SKSS-MD2003)



Spirometer Various Manufacturer by LABYCARE

Description :

Spirometer is a medical device used to measure lung function, including volume and airflow. It helps diagnose respiratory conditions such as asthma and COPD. Available in manual and digital models, it is widely used in hospitals, clinics, and laboratories for pulmonary function testing and patient monitoring.

Specifications :

Available in manual (water-sealed) and digital models
Measures lung volumes: tidal volume, vital capacity, FVC, FEV1
Clear display (digital) or calibrated scale (manual)
High accuracy and repeatability
User-friendly interface with simple operation

Flow/volume measurement with sensors

labycare

www.labycare.com | Email: labycare@gmail.com