

Labycare-Operation Table Various (SKSS-MD2005)



Operation Table Various Manufacturer by LABYCARE

Description :

Operation table is a versatile surgical platform designed to support patients during various procedures. It offers adjustable positioning, stability, and comfort with advanced control systems. Suitable for multiple surgeries, it ensures precise positioning, durability, and ease of operation, making it essential equipment in operation theatres and surgical centers.

Specifications :

- Table Type: Hydraulic / Electric
- Height Adjustment: Approx. 700 – 1000 mm
- Table Top: Radiolucent, sectional design
- Positioning: Trendelenburg / Reverse / Lateral tilt
- Load Capacity: ? 150 – 200 kg
- Base & Frame: Stainless steel, heavy-duty construction

labycare

www.labycare.com | Email: labycare@gmail.com