


Labycare-Perfusion Assembly (SKSS-MD2006)



Perfusion Assembly Manufacturer by LABYCARE

Description :

Perfusion assembly is a laboratory apparatus used to study the physiological effects of drugs on isolated tissues or organs. It ensures continuous supply of nutrient solution under controlled conditions. Commonly used in pharmacology experiments, it helps maintain tissue viability and supports accurate observation of biological responses.

Specifications :

Material: Glass / Stainless Steel
Flow Control: Adjustable
Usage: Tissue/organ perfusion
Compatibility: Organ bath system
Mounting: Stand compatible

Documents : <images/bbbcb326b491f4f858a6ee61f49eb3aa.png>

labycare

www.labycare.com | Email: labycare@gmail.com