

## Labycare-Heart Levers Various (SKSS-MD2007)



Heart Levers Various Manufacturer by LABYCARE

### **Description :**

Heart levers are laboratory instruments used in physiology experiments to record heart movements of animals on a kymograph. They provide accurate mechanical amplification of cardiac contractions for study and analysis. Commonly used in educational and research laboratories, they help demonstrate cardiovascular physiology effectively.

### **Specifications :**

Material: Brass / Stainless Steel

Type: Various (Starling / Simple)

Mounting: Stand compatible

Application: Cardiac physiology experiments

Finish: Polished, corrosion-resistant

**Documents :** [images/7a1bd2b61d3a7233e12c511d7b4920e2.png](#)

---

labycare

www.labycare.com | Email: labycare@gmail.com