

# Labycare-Brodie Universal Levers (SKSS-MD2008)



Brodie Universal Levers Manufacturer by LABYCARE

## **Description :**

Brodie universal lever is a precision instrument used in physiology experiments to record movements of muscles, heart, or organs. It amplifies small mechanical changes into measurable recordings on kymographs. Widely used in laboratories, it ensures accurate and sensitive recording of biological responses under controlled experimental conditions.

## **Specifications :**

Material: Aluminium / Stainless Steel

Type: Universal Lever

Mounting: Stand compatible

Usage: Muscle/organ movement recording

Finish: Polished and rust-resistant

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)