

Labycare-Pulleys Various (SKSS-MD2009)



Pulleys Various Manufacturer by LABYCARE

Description :

Laboratory pulleys are used in physiology and physics experiments to guide threads or wires smoothly while recording movements. They help reduce friction and ensure accurate transmission of motion in experimental setups. Commonly used with kymographs and levers, they support precise and controlled observation of mechanical responses.

Specifications :

Material: Metal / Plastic

Type: Single / Double pulley

Mounting: Clamp/stand compatible

Rotation: Free-moving

Use: Motion transmission in experiments

labycare

www.labycare.com | Email: labycare@gmail.com