

## Labycare-Baths Various ( SKSS-MD2010)



Baths Various Manufacturer by LABYCARE

### **Description :**

Medical baths are devices used to maintain controlled temperature conditions for samples, instruments, or therapeutic applications. Commonly used in laboratories and healthcare settings, they ensure uniform heating and precise temperature regulation. Designed for reliability and safety, they support various medical, diagnostic, and research procedures efficiently.

### **Specifications :**

Temperature Range: Ambient to 100°C (or as required)

Temperature Control: Digital / Thermostatic

Capacity: Various (2L to 20L or more)

Material: Stainless steel inner chamber

Heating Type: Electric immersion heater

Display: LED / Digital temperature display

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)