

Labycare-Piston Recorder Various (SKSS-MD2012)



Piston Recorder Various Manufacturer by LABYCARE

Description :

Piston recorder is a medical laboratory device used to record physiological movements such as respiration or pressure changes. It converts mechanical motion into graphical recordings for analysis. Commonly used in physiology experiments, it ensures accurate and continuous recording of biological responses in controlled laboratory conditions.

Specifications :

Type: Piston-operated recording device

Material: Glass cylinder with metal fittings

Movement: Smooth, frictionless linear piston motion

Recording: Compatible with kymograph/drum recording system

Sensitivity: High sensitivity for small pressure changes

Mounting: Stand compatible with clamp support

Application: Respiration/pressure recording in physiology experiments

Documents : <images/57199b21a6e9eb18ca25526fdc776172.png>

labycare

www.labycare.com | Email: labycare@gmail.com