

## Labycare-Dubios Reymond Key (SKSS-MD2017)



Dubios Reymond Key Manufacturer by LABYCARE

### **Description :**

Du Bois-Reymond key is a laboratory switching device used in physiology experiments to control electrical stimulation circuits. It allows making or breaking the circuit instantly, ensuring precise delivery of stimuli to tissues. Commonly used with induction coils, it supports accurate and controlled experimental observations.

### **Specifications :**

Type: Electrical contact key (make-and-break switch)

Material: Brass contacts with ebonite/insulated base

Operation: Manual press key for circuit make/break

Contact Type: Spring-loaded, quick response

Mounting: Base mounted / bench type

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)