

Labycare-Press Contact Key (SKSS-MD2018)



Press Contact Key Manufacturer by LABYCARE

Description :

Press contact key is a simple electrical switch used in physiology experiments to make or break a circuit instantly. It enables precise control of electrical stimulation in laboratory setups. Commonly used with stimulators and induction coils, it ensures accurate timing and reliable performance during experimental procedures.

Specifications :

Type: Press contact switch

Material: Brass contacts, insulated base

Operation: Manual press (momentary contact)

Terminals: Binding posts

Mounting: Bench/base mounted

Documents : <images/92f78d5b59ed9bd5f5f389f4976932f9.png>

labycare

www.labycare.com | Email: labycare@gmail.com