

Labycare-Simple Electrode (SKSS-MD2019)



Simple Electrode Manufacturer by LABYCARE

Description :

Simple electrode is a basic conducting device used in physiology and laboratory experiments to deliver or receive electrical signals from tissues. It is commonly used for stimulation or recording purposes. Designed for easy handling, it ensures stable contact and reliable performance in experimental and diagnostic applications.

Specifications :

Material: Copper / Silver / Stainless Steel

Type: Surface / Needle electrode

Connection: Wire with terminal clips

Usage: Stimulation/recording

Finish: Corrosion-resistant

Reusable: Yes

labycare

www.labycare.com | Email: labycare@gmail.com