

## Labycare-Tuning Fork N-100 (SKSS-MD2020)



Tuning Fork N-100 Manufacturer by LABYCARE

### **Description :**

Tuning Fork N-100 is a medical diagnostic instrument used to assess hearing and vibration sense. It produces a fixed frequency sound when struck, aiding in neurological and auditory examinations. Commonly used in clinics and hospitals, it ensures accurate and reliable evaluation of patient sensory responses.

### **Specifications :**

Type: Medical tuning fork

Frequency: 100 Hz

Material: Aluminium alloy / Steel

Finish: Polished, corrosion-resistant

Use: Hearing and vibration testing

Reusable: Yes

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)