

Labycare-Hook & Weight Set 11gm (SKSS-MD2021)



Hook & Weight Set 11gm Manufacturer by LABYCARE

Description :

Hook and weight set (11 gm) is used in physiology experiments to apply controlled load on tissues or levers. It helps measure muscle contraction and responses accurately. Designed for precision and ease of use, it is commonly used in laboratory setups with kymographs and recording systems.

Specifications :

Weight: 11 gm

Material: Brass / Stainless Steel

Type: Hook with calibrated weight

Use: Muscle/lever loading

Finish: Polished, rust-resistant

Reusable: Yes

labycare

www.labycare.com | Email: labycare@gmail.com