

Labycare-Hook & Weight Set 110gm (SKSS-MD2022)



Hook & Weight Set 110gm Manufacturer by LABYCARE

Description :

Hook and weight set (110 gm) is used in physiology experiments to apply controlled load on tissues, muscles, or levers. It assists in studying contraction and response under different loads. Designed for accuracy and durability, it is suitable for laboratory setups involving kymographs and recording systems.

Specifications :

Weight: 110 gm

Material: Brass / Stainless Steel

Type: Hook with weight set

Use: Muscle/lever loading

Finish: Polished, rust-resistant

Reusable: Yes

labycare

www.labycare.com | Email: labycare@gmail.com