

Labycare-Muscle Grip (SKSS-MD2023)



Muscle Grip Manufacturer by LABYCARE

Description :

Muscle grip is a laboratory instrument used in physiology experiments to hold muscle tissues securely during recording or stimulation. It ensures firm attachment without damaging the tissue, allowing accurate observation of contractions. Commonly used with kymographs and levers, it supports reliable and controlled experimental procedures.

Specifications :

Material: Stainless Steel / Brass

Type: Tissue holding clamp

Grip: Serrated / smooth jaws

Mounting: Stand compatible

Use: Muscle fixation in experiments

Finish: Polished, rust-resistant

labycare

www.labycare.com | Email: labycare@gmail.com