

## Labycare-Pithing Needle S. S (SKSS-MD2024)



Pithing Needle S. S Manufacturer by LABYCARE

### **Description :**

Pithing needle is a stainless steel instrument used in physiology and biological experiments for humane destruction of brain and spinal cord in laboratory animals. It ensures quick and effective pithing. Designed with a sharp tip and sturdy handle, it provides precision, control, and safety during procedures.

### **Specifications :**

Material: Stainless Steel

Type: Pithing/Dissection needle

Tip: Sharp, pointed

Handle: Metal / insulated grip

Use: Physiology/dissection procedures

Finish: Polished, rust-resistant

---

labycare

www.labycare.com | Email: labycare@gmail.com