

Labycare-Fixing Pins (SKSS-MD2025)



Fixing Pins Manufacturer by LABYCARE

Description :

Connecting wires are insulated electrical leads used in laboratory and medical experiments to connect instruments and circuits. They ensure safe and efficient transmission of electrical signals. Commonly used with stimulators, power supplies, and recording devices, they provide reliable connectivity and flexibility in experimental setups.

Specifications :

Type: Insulated connecting leads

Conductor: Copper wire

Insulation: PVC / Silicone

Length: Approx. 1–1.5 meters

Connectors: Banana plug / Crocodile clips

Color: Red/Black (standard)

Voltage Rating: Low voltage lab use

labycare

www.labycare.com | Email: labycare@gmail.com