

## Labycare-Muscle Chamber (SKSS-MD2026)



Muscle Chamber Manufacturer by LABYCARE

### **Description :**

Muscle chamber is a laboratory apparatus used to hold isolated muscle or tissue in physiological solution during experiments. It maintains proper conditions for tissue viability and allows recording of contractions. Commonly used in pharmacology and physiology labs, it supports accurate study of muscle responses and drug effects.

### **Specifications :**

Material: Glass / Acrylic / Stainless Steel  
Capacity: Approx. 20–50 mL  
Design: With inlet/outlet for solution flow  
Mounting: Stand compatible  
Accessories: Tissue holder, aeration tube

Use: Isolated tissue/muscle experiments

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)