

Labycare-Mygraph Board With Crank Lever (SKSS-MD2027)



Mygraph Board With Crank Lever Manufacturer by LABYCARE

Description :

Mygraph board with crank lever is used in physiology experiments to study muscle contraction and movements. It provides a stable platform for mounting tissue and recording devices. The crank lever mechanism helps apply controlled motion, enabling accurate recording of muscle responses in laboratory experiments.

Specifications :

Material: Wooden/metal base

Mechanism: Crank lever system

Mounting: Stand and clamp compatible

Operation: Manual rotation

Use: Muscle contraction experiments

Finish: Durable and corrosion-resistant parts

labycare

www.labycare.com | Email: labycare@gmail.com