


Labycare-White Glazed Paper Various (SKSS-MD2028)



White Glazed Paper Various Manufacturer by LABYCARE

Description :

White glazed paper is used in laboratory experiments for recording physiological tracings and observations. It provides a smooth, glossy surface suitable for clear markings with ink or stylus. Commonly used with kymographs and recording devices, it ensures accurate and visible graphical representation of experimental data.

Specifications :

Type: Glazed recording paper

Surface: Smooth, glossy finish

Color: White

Size: Various (sheets/rolls)

Compatibility: Kymograph/recording systems

Use: Graphical recording and documentation

labycare

www.labycare.com | Email: labycare@gmail.com