

Labycare-Perimeter Lister (SKSS-MD2029)



Perimeter Lister Manufacturer by LABYCARE

Description :

Perimeter Lister is an ophthalmic instrument used to measure the visual field of the eye. It helps detect defects in peripheral vision caused by glaucoma or neurological disorders. Designed for accuracy and ease of use, it is commonly used in eye clinics and hospitals for diagnostic evaluation.

Specifications :

Type: Manual visual field perimeter

Arc: Graduated semicircular arc

Material: Metal with matte finish

Measurement: Degrees of visual field

Mounting: Table/stand mounted

Use: Visual field testing in ophthalmology

labycare

www.labycare.com | Email: labycare@gmail.com