

Labycare-Out Flow Recorders Various (SKSS-MD2030)



Out Flow Recorders Various Manufacturer by LABYCARE

Description :

Outflow recorder is a laboratory instrument used to measure and record the rate of fluid flow from an organ or system during experiments. It helps study physiological responses and drug effects by providing accurate flow measurements, commonly used in pharmacology and physiology laboratory setups.

Specifications :

Type: Fluid outflow recording apparatus
Material: Glass chamber with metal fittings
Measurement: Continuous fluid flow rate recording
Scale: Graduated for volume measurement
Connection: Inlet/outlet tubing ports
Compatibility: Kymograph/recording system
Mounting: Stand/clamp support

Accuracy: High sensitivity for small flow changes

labycare

www.labycare.com | Email: labycare@gmail.com