


Labycare-Commutators Various (Product Cat # SKSS-MD2031)



Commutators Various Manufacturer by LABYCARE

Description :

Commutators are laboratory switching devices used to reverse or control the direction of electrical current in experimental circuits. Commonly used in physiology setups, they help manage stimulation and recording connections. Designed for precision and reliability, they ensure smooth switching and accurate control during laboratory experiments.

Specifications :

Type: Rotary / slide commutator

Material: Brass contacts with insulated base

Operation: Manual switching

Terminals: Multiple binding posts

Function: Current direction control

Mounting: Base/bench type

Compatibility: Lab circuits and stimulators

labycare

www.labycare.com | Email: labycare@gmail.com