

Labycare-Bycycle (Martin) Ergograph (SKSS-MD2033)



Bicycle (Martin) Ergograph Manufacturer by LABYCARE

Description :

Bicycle (Martin) ergograph is a physiological instrument used to measure work done and fatigue of muscles during exercise. It simulates cycling activity under controlled conditions, allowing study of muscular performance. Commonly used in physiology laboratories, it helps analyze endurance, workload, and fatigue patterns accurately.

Specifications :

Type: Mechanical ergograph (cycle type)

Measurement: Work output and fatigue

Mechanism: Pedal-driven resistance system

Adjustment: Variable load control

Display: Mechanical/dial reading

Mounting: Floor/bench type

Documents : [images/1a6710fbaa7351a5019953fbd7667049.png](#)

labycare

www.labycare.com | Email: labycare@gmail.com