

Labycare-Mosses Ergograph (SKSS-MD2034)



Mosses Ergograph Manufacturer by LABYCARE

Description :

Mosso's ergograph is a physiology instrument used to measure muscle work and fatigue, typically using finger contractions. It records the amount of work performed over time, helping study endurance and fatigue. Widely used in laboratories, it provides accurate analysis of muscular performance under controlled experimental conditions.

Specifications :

Type: Mechanical ergograph

Measurement: Muscle work and fatigue

Mechanism: Finger contraction with weight load

Recording: Compatible with kymograph

Adjustment: Variable load system

Mounting: Table/stand mounted

labycare

www.labycare.com | Email: labycare@gmail.com