

Labycare- Newton Colour Disc (SKSS-MD2035)



Newton Colour Disc Manufacturer by LABYCARE

Description :

Newton colour disc is a physics and physiology demonstration device used to show that white light is composed of different colors. When rotated rapidly, the colored segments blend to appear white. It is commonly used in laboratories and classrooms to demonstrate principles of color mixing and light.

Specifications :

Type: Rotating color disc

Material: Plastic / metal disc with color segments

Operation: Manual / motorized rotation

Mounting: Stand/axle mounted

Use: Color mixing and optics demonstration

labycare

www.labycare.com | Email: labycare@gmail.com