

Labycare-Electrical Colour Mixer (SKSS-MD2037)



Electrical Colour Mixer Manufacturer by LABYCARE

Description :

Electrical colour mixer is a motorized device used to demonstrate color mixing by rotating colored discs at controlled speeds. It enables uniform blending of colors to produce white or intermediate shades. Commonly used in laboratories and classrooms, it provides accurate and repeatable results for studying light and color principles.

Specifications :

Type: Motorized color mixer

Operation: Electric motor driven

Speed Control: Adjustable rotation speed

Material: Metal/plastic disc with color segments

Power Supply: 220–230V AC, 50/60 Hz

Mounting: Table/stand mounted

labycare

www.labycare.com | Email: labycare@gmail.com