

Labycare-Large Metal Screen (SKSS-MD2038)



Large Metal Screen Large Metal Screen

Description :

Large metal screen is used in laboratory experiments for projecting and observing light patterns, shadows, and images. It provides a flat, stable surface with matte finish for clear visibility. Commonly used in optics and physics labs, it ensures accurate observation and supports demonstration of light and projection principles effectively.

Specifications :

Material: Mild steel / Aluminium

Surface Finish: Matte white / non-reflective coating

Size: Large (various dimensions available)

Mounting: Stand-mounted / wall mount option

Frame: Rigid metal frame for stability

Portability: Lightweight, easy to move

labycare

www.labycare.com | Email: labycare@gmail.com