

Labycare-Hair Aesthesiometer (Von Freys) (SKSS-MD2039)



Hair Aesthesiometer (Von Freys) Manufacturer by LABYCARE

Description :

Hair aesthesiometer (Von Frey's) is a diagnostic instrument used to assess touch and pressure sensitivity of the skin. It consists of calibrated monofilaments that apply specific forces. Commonly used in neurological and physiological studies, it helps evaluate sensory nerve function and detect abnormalities accurately.

Specifications :

Type: Monofilament sensory testing device
Material: Nylon filaments with handles
Range: Calibrated force levels (various)
Measurement: Touch/pressure sensitivity
Set: Multiple filaments of different thickness

Use: Neurological and sensory testing

labycare

www.labycare.com | Email: labycare@gmail.com