

Labycare-Simple Hair Aesthesiometer (SKSS-MD2040)



Simple Hair Aesthesiometer Manufacturer by LABYCARE

Description :

Simple hair aesthesiometer is a basic instrument used to assess touch sensitivity of the skin. It uses a single calibrated filament to apply a controlled stimulus. Commonly used in physiology and clinical settings, it helps evaluate sensory nerve response in a simple and effective manner.

Specifications :

Type: Single filament aesthesiometer

Material: Nylon hair with handle

Force: Fixed calibrated pressure

Use: Touch sensitivity testing

Design: Handheld, lightweight

Reusable: Yes

labycare

www.labycare.com | Email: labycare@gmail.com