

Labycare-Heat & Cold Spot Hammer (SKSS-MD2042)



Heat & Cold Spot Hammer Manufacturer by LABYCARE

Description :

Heat and cold spot hammer is a diagnostic instrument used to test thermal sensitivity of the skin. It has two ends designed for hot and cold stimulation. Commonly used in physiology and neurological examinations, it helps identify temperature sensation and detect sensory nerve abnormalities accurately.

Specifications :

Type: Thermal sensory testing instrument

Material: Metal with insulated handle

Ends: Dual (heat and cold)

Use: Temperature sensation testing

Design: Handheld, lightweight

Application: Physiology and neurological exams

labycare

www.labycare.com | Email: labycare@gmail.com